

SK-L407

5W WIRELESS HEADLIGHT




SK-L407

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Voltage	AC 90-240V 50-60Hz
▪ Blub Power	LED 5W
▪ Color temperature	4000~5000K
▪ Bulb life	≥50000h
▪ Illumination	≥30000Lux
▪ DC charging time	≤5h
▪ Power supply time	≥4h+4h
▪ Load weight on head	330g
▪ Net weight	480g
▪ Package size	230x230x140mm
▪ Package material	bandbox



Description

It adopt 5W LED as the light source, with high brightness and CRI, long life span.

Light distribution uniform, light spot clear and bright. Wireless, two batteries, easy and quickly to replace, work is guaranteed. There is a joint under the head, could attach loupe, more suitable for difficult operation.

Two batteries, light spot size and brightness adjustable.

Application

Surgery (cardiac surgery, neurosurgery, brain surgery), ENT, Dental, Plastic surgery, Pet medical.

Features

LED source.
High CRI, high illumination.
Brightness adjustable.



Magnifying glass

It can be used with magnifying glass.



Light spot adjustable



The size of the spot is adjustable



Dual battery supply

Task guaranteed.



Optional

Optional aluminum alloy packaging.
High-end and classy.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🕒 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue