



A045-1

GYNECOLOGICAL TABLE

Registration code: SKF-A





SKSurgitech[™]
by Saikang Medical

Care you precisely

A045-1

FEATURES AND HIGHLIGHTS



Technical parameters

| | |
|-----------------------|-----------------|
| External size (LxWxH) | 1820x500x800 mm |
| Back-rest adjustment | 0-55° (±10°) |
| Knee-rest adjustment | 0-85° (±10°) |

Technical configuration

| | |
|--|-------|
| <input checked="" type="checkbox"/> Lithotomy Rod | 1set |
| <input checked="" type="checkbox"/> Hand Bar | 1set |
| <input checked="" type="checkbox"/> Anti-slip Bed Feet Cover | 4pcs |
| <input checked="" type="checkbox"/> Mechanical Level | 2sets |
| <input checked="" type="checkbox"/> Stainless Steel Basin | 1pc |



Design

Simple structure, practical; Whole stainless steel material meets the needs of the surgical room. This bed is designed for gynecological diagnosis, operation and parturition.



Material

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth wedding.



Leather

Flame retardant, yellow resistance, mildew resistance, bending resistance, 100,000 times wear resistance.



Sponge

Flame retardant, 4 cm thickness, high-density foam, density 30 kg/m³.



Lithotomy Rod

Lithotomy Rod, adjustable in height and different angles.



Hand Bar

Granular surface to increase friction and prevent hand slip.



Stainless Steel Basin

Removable stainless steel basin, convenient to take out.



Feet Cover

4 pcs Anti-slip bed feet covers.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🕒 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue