



A04

ORTHOPEDIC TRACTION FRAME

Registration code: SKT-C



SKSurgitech[™]
by Saikang Medical

Care you precisely

A04

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Traction length	2000mm
▪ Traction width	600mm
▪ Allowed can equipped table top minimum height	≤740
▪ Movements	Model:1005
▪ Footshoes displacement	≥200mm
▪ Leg support minimum height	≤540mm
▪ Leg support elevation range	≥200mm
▪ Adjustment of the foot shoes' angle range	≥60°
▪ The angle range of main horizontal support frame	≥180°

Technical configuration

▪ Lift crank handle	1piece
▪ Lateral support	1piece
▪ Knee support	1piece
▪ Straight cvlinderholder	1piece





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue