



A100-4A

ELECTRIC OPERATING TABLE

Registration code: SKL-A



SKSurgitech[™]
by Saikang Medical


Care you precisely

A100-4A

FEATURES AND HIGHLIGHTS

Technical parameters

External Size (LxWxH):	2060x500x710-1010 mm
Built-in Kidney Bridge:	0-110 mm
Back Plate Up/Down:	80°/20°
Head Plate Up/Down:	40°/90°
Leg Plate Up/Down:	15°/90°
Trendelenburg/Reverse Trendelenburg:	25°/25°
Lateral Tilt Left/Right:	15°/15°
Horizontal Sliding:	350 mm
Battery Standby Time:	80 h
Power Supply:	220V 50HZ 350VA
Power Control Circuit:	12V
Plug:	as required
X-rays translucent:	YES



- C-arm compatible: YES
- Multi-functional orthopedic frame: YES

Technical configuration

<input checked="" type="checkbox"/> Bed Frame	1pc
<input checked="" type="checkbox"/> Anesthesia Screen Frame	1pc
<input checked="" type="checkbox"/> Mattress	1set
<input checked="" type="checkbox"/> Kidney Bridge Handle	1pc
<input checked="" type="checkbox"/> Arm Support with Cushion	2pcs
<input checked="" type="checkbox"/> Leg Support with Cushion	2pcs
<input checked="" type="checkbox"/> Shoulder Support with Cushion	2pcs
<input checked="" type="checkbox"/> Body Support with Cushion	2pcs
<input checked="" type="checkbox"/> Foot Support with Cushion	2pcs
<input checked="" type="checkbox"/> Hand Controller	1pc

Remote Panel & Column Panel

Leg Plate With Imported Gas Spring For Easy Control.



Kidney Bridge

Built-in kidney bridge, more convenient and nice appearance.



Wheel

6" castors import from germany.



Leg Plate

Leg Plate With Imported Gas Spring For Easy Control.



Material

The base seat and upright column covering are made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned.



Braking System

Pedal Brake, Easy To Move And Lock.



Powder Supply

Internal battery maintaining power supply for 50-80 operations and with low battery warning function.



Shoulder Support with Cushion



Body Support with Cushion



Arm Support with Cushion



Leg Support with Cushion



Anesthesia Screen Frame



Lowest Position of Table Top
710mm



Highest Position of Table Top
1010mm




Back Plate Up
0°-80°



Back Plate Down
0°-20°



Reversed Trendelenburg
0°-25°



Trendelenburg
0°-25°



Lateral Tilt Left
0°-15°



Lateral Tilt Right
0°-15°



Horizontal Sliding
350mm



Head Plate Up
0°-40°



Head Plate Down
0°-90°



Leg Plate Rest
-90°~15°





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa





Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 ·  24/7
 export@sksurgitech.com
 www.sksurgitech.com



Company video



General catalogue