



A100-4

ELECTRIC OPERATING TABLE

Registration code: SKL-A




SKSurgitech[™]
by Saikang Medical

Care you precisely

A100-4

FEATURES AND HIGHLIGHTS



Technical parameters

External Size (LxWxH):	2070x550x700-1000 mm
Built-in Kidney Bridge:	0-120 mm
Back Plate Up/Down:	75°/20°
Head Plate Up/Down:	45°/90°
Leg Plate Up/Down:	15°/90°
Trendelenburg/Reverse Trendelenburg:	25°/25°
Lateral Tilt Left/Right:	15°/15°
Horizontal Sliding:	300 mm
X-rays translucent:	YES
C-arm compatible:	YES
Multi-functional orthopedic frame:	YES

Technical configuration

<input checked="" type="checkbox"/> Bed Frame	1pc
<input checked="" type="checkbox"/> Anesthesia Screen Frame	1pc
<input checked="" type="checkbox"/> Mattress	1set
<input checked="" type="checkbox"/> Kidney Bridge Handle	1pc
<input checked="" type="checkbox"/> Arm Support with Cushion	2pcs
<input checked="" type="checkbox"/> Leg Support with Cushion	2pcs
<input checked="" type="checkbox"/> Shoulder Support with Cushion	2pcs
<input checked="" type="checkbox"/> Body Support with Cushion	2pcs
<input checked="" type="checkbox"/> Foot Support with Cushion	2pcs

Controller

All function controlled by hand controller, and can be locked.



C-arm

Extra-long horizontal sliding, available for c-arm operation.




Kidney Bridge

Built-in kidney bridge, more convenient and nice appearance.



Wash Basin cutout

Five section table top including Head, Back, Shoulder, Hip and Leg, all section applicable for X ray use and under each section with X ray cassette layer.



Head Plate

Detachable double head plate and leg plate meet operation position needs and use Germany cylinder assist system.



Power

Internal battery maintaining power supply for 50-80 operations and with low battery warning function.



Shoulder Support with Cushion



Body Support with Cushion



Arm Support with Cushion



Leg Support with Cushion



Lowest Position of Table Top
700mm



Highest Position of Table Top
1000mm



Longitudinal Adjustment
0-300mm



Longitudinal Adjustment
0-300mm



Reversed Trendelenburg
0°-22°



Trendelenburg
0°-25°



Lateral Tilt Left
0°-20°



Lateral Tilt Right
0°-20°



Kidney Bridge Elevation
110mm



Anorectum Operation
Position



E.N.T Operation Position



(Optional) Head Rest



(Optional) Lateral Arm
Support



(Optional) Waste Basin For
Gynecology



(Optional) Orthopedics Traction Frame






CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🕒 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue