



A106-2

ELECTRIC-HYDRAULIC OPERATING TABLE

Registration code: SKT-B



SKSurgitech[™]
by Saikang Medical

Care you precisely

A106-2

FEATURES AND HIGHLIGHTS



Technical parameters

External Length:	2080mm
External Width:	550mm
Table Height Min:	500mm
Table Height Max:	800mm
Reverse Trendelenburg:	0°-25°
Trendelenburg:	0°-25°
Lateral Tilt Left/Right:	0°-20°/0°-20°
Head Plate Up/Down:	0°-45°/0°-90°
Back Plate Up/Down:	0°-80°/0°-20°
Foot Plate Up/Down:	0°-15°/0°-90°

Foot Plate Horizontal Positioning:	0°-90°
Patient Weight Capacity:	250kg
X-rays translucent:	YES
C-arm compatible:	NO
Multi-functional orthopedic frame:	YES

Technical configuration

<input checked="" type="checkbox"/> Anaesthesia Screen	1pc
<input checked="" type="checkbox"/> Hand Remote	1pc
<input checked="" type="checkbox"/> Kidney Bridge Handle	1pc
<input checked="" type="checkbox"/> Arm Support	1pc
<input checked="" type="checkbox"/> Foot Support	1pc
<input checked="" type="checkbox"/> Mattress	1set
<input checked="" type="checkbox"/> Shoulder Support	1pair
<input checked="" type="checkbox"/> Long Fixing Clamp	1pair
<input checked="" type="checkbox"/> Fixing Clamp	7pcs
<input checked="" type="checkbox"/> Power Line	1pc



Material

The base, column and side hoder are all made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned.



Powder Supply

Internal battery maintaining power supply for 50-80 operations and with low battery warning function.




Remote Panel

All function controlled by hand controller in the column, and can be locked.



Hydraulic System

Brand hydraulic system (German Hoerbiger) make table with low noise, high safety and reliability.




Control System

Foot treadle brake device which is easily for moving operating table.



Super Low Position

Super low position (500mm) , suitable for operation on head.



Kidney Bridge

Built-In Kidney Bridge, More Convenient And Nice Appearance.



Screen



Armrest



Knee Curtch



Shoulder Support



Body Support





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa





Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 ·  24/7
 export@sksurgitech.com
 www.sksurgitech.com



Company video



General catalogue