



# A2000

## ELECTRIC OPERATING TABLE

Registration code: SKL-A



**SKSurgitech**<sup>TM</sup>  
by Saikang Medical

Care you precisely

# A2000

## FEATURES AND HIGHLIGHTS



### Technical parameters


External size (LxWxH):	2100x480x(680-930)mm
Reverse trendelenburg:	0-25°
Trendelenburg:	0-25°
Lateral tilt left/right:	0-15°
Head plate up/down:	45°/90°
Back plate up/down:	70°/15°
Leg plate up/down:	45°/90°
Leg plate outwards:	0-90°
Kidney bridge elevation:	0-100mm
Horizontal sliding:	0-350mm
Safe working load:	150kg

### Technical configuration

<input checked="" type="checkbox"/> Bed frame	1 pc
<input checked="" type="checkbox"/> Mattress	1 set
<input checked="" type="checkbox"/> Anesthesia screen frame	1 pc
<input checked="" type="checkbox"/> Shoulder holder	2 pcs
<input checked="" type="checkbox"/> Arm holder	2 pcs
<input checked="" type="checkbox"/> Leg support	2 pcs
<input checked="" type="checkbox"/> Foot support	2 pcs
<input checked="" type="checkbox"/> Hand controller	1 pc
<input checked="" type="checkbox"/> Sliding block	10 pcs
<input checked="" type="checkbox"/> Spanner	4 pcs

**Material**

The base material for the table is high thickness 304 stainless steel, seamless welding by robot, with modern appearance and more stable.



**Mattress**

Removable memory mattress, good rebound resilience, anti-static electricity and anti-microbial.



**Brake**

Pedal brake, easy to move and lock.



**Safety**

The table is provided with an emergency stop switch on the base column, help to stop the operating table running in case of emergency situation, ensure the safety of operation.



**Leg plate**

Leg plate with imported gas spring for easy control.



**Kidney bridge**

Manual kidney bridge, more convenience.



**Battery could support**

Battery could support about 50-80 complete operation cycles. Working time could be last for one week according to useage and help to keep the operating table running in the event of power failure.



**Armrest**



**Anesthesia Frame**



**Leg Crutches**



**Shoulder Support**



**Body Supporter**



Professional functions

Lowest position of table top

**680**mm



Highest position of table top

**930**mm



Back plate up

**0°-70°**



Back plate down

**0°-15°**



Reversed trendelenburg

**0°-25°**



Trendelenburg

**0°-25°**



Right lateral tilt

**0°-15°**



Left lateral tilt

**0°-15°**



Head plate upward

**0°-45°**



Head plate downward

**0°-90°**



Leg section evagination

**0°-90°**



Leg plate up/down

**-45°-90°**



Kidney bridge elevation

**0-100**mm





## **CONSTANT INNOVATION**



Design



Research



Manufacturing



Maintenance



Service

## **WORLDWIDE BRANCHES**

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa




**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
 +86 18021231901 ·  24/7  
 export@sksurgitech.com  
 www.sksurgitech.com



Company video



General catalogue