

A98-3


ELECTRIC OBSTETRIC BED

Registration code: SKL-A



A98-3

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External size (LxWxH)	2140 x 1120 x (595-905) ±10mm
▪ PP back section	775x875mm
▪ PP seat section	430x880mm
▪ PP leg section	650x794mm
▪ Back-rest adjustment	0-70°
▪ Trendelenburg	0-10°


Technical configuration

<input checked="" type="checkbox"/> Leg Holder (PU material)	1pair
<input checked="" type="checkbox"/> Grab Handle	1pair
<input checked="" type="checkbox"/> Basin	1pc
<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> IV Pole Prevision	2pcs
<input checked="" type="checkbox"/> Controller	1pc
<input checked="" type="checkbox"/> Satinless Steel Guardrail	1pair
<input checked="" type="checkbox"/> 6" Castor	4pcs
<input checked="" type="checkbox"/> 80mm Mattress	1set




Material

The base seat and height adjustment system are made of powder coated steel, high gloss, anti corrosion, easy to clean.



Bed board

Bed surface has PP cover, made by virgin material. The PP cover is detachable, and easy to clean.



Head&Foot Board

The bed platform and head&foot board made by virgin PP material, the head&foot board is movable.



Hand controller

Equipped with one hand controller.




Guardrail

Handrail is rotatable, when you don't use it, you can put down the guardrail.




Handle

The hand bar can be dismantled and movable.




Lithotomy rod

The leg support can be dismantled and movable.



Foot support

The foot support can be dismantled and movable.



Castor

6 inch central controlled castor for moving obstetric table.



Accessory

Equipped with stainless steel sewage basin and other various accessories.

Mattress

The mattress adopts high density memory sponge and is once formed. Mattress is movable.

Foot section

foot section can be retracted into the bed, convenient for gynecological and delivery surgical.



External Size (LxW)

2140x1120 mm



Back-rest adjustment

70°



Height adjustment

595-905 mm



Trendelenburg

10°



Foot section into the bed



CONSTANT INNOVATION

Design

Research

Manufacturing

Maintenance

Service



WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd

Sala Thammasop Subdistrict
Thawi
Watthana District
Bangkok Thailand
TEL: +66 642741046

UK Branch

Uk Healthward International Limited

DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +31 (0) 2021111106

Russia Branch

LLC 'SAIKANG'

454004, Russian Federation,
Chelyabinsk,
Universitetskaya Naberezhnaya,
Russia
TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC

Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

 +86 18021231901 ·  24/7

 export@sksurgitech.com

 www.sksurgitech.com



Company video



General catalogue