


A98-4

ELECTRIC OBSTETRIC BED



A98-4

FEATURES AND HIGHLIGHTS




Technical parameters

▪ External Size (LxW)	1980x960mm(±30mm)
▪ Height	650-940mm(±20mm)
▪ Back-rest Adjustment	0-60°
▪ Reverse Trendelenburg	0-7°
▪ Leg-rest Adjustment	-9°-35°, leg section is removable
▪ Weight capacity	200KG



Technical configuration

<input checked="" type="checkbox"/> Linear motor	3pcs
<input checked="" type="checkbox"/> 5" Double side castors, central locking	4pcs
<input checked="" type="checkbox"/> Central locking pedal	2pcs
<input checked="" type="checkbox"/> ABS head board with safe lock	1pc
<input checked="" type="checkbox"/> ABS side rail with angle indicator	1set
<input checked="" type="checkbox"/> Hand grip	1pair
<input checked="" type="checkbox"/> Leg holder	1pair
<input checked="" type="checkbox"/> Basin	1pc
<input checked="" type="checkbox"/> Hand controller	1pc
<input checked="" type="checkbox"/> Power line	1pc
<input checked="" type="checkbox"/> Leg strip	1 pair
<input checked="" type="checkbox"/> IV pole prevision	2pcs



Material

The base seat and height adjustment system are made of powder coated steel, high gloss, anti corrosion, easy to clean.



Bed surface

Bed surface by plastic steel, made by virgin material. The cover is detachable, and easy to clean.



Head Board and Guardrail

The Head board and side rail made by virgin ABS material. The head board is movable.



Hand controller

Equipped with one hand controller.




Handle

The hand bar can be dismantled and movable.



Lithotomy rod

The leg support can be dismantled and movable.



Caster

5 inch central controlled caster for moving obstetric table. TPR tire, quiet and wear-resistant, can be washed and sterilized by high temperature, save anti-winding hard shell, united forming without bolts.



Accessory

Equipped with stainless steel sewage basin and other various accessories.



Mattress

The mattress adopts high density memory sponge and is once formed. Mattress is movable.



Foot section

Foot section can be retracted, convenient for gynecological and delivery surgical.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue