



A99-7

ELECTRIC  
OBSTETRIC BED

Registration code: SKL-A





SKSurgitech™  
by Saikang Medical

Care you precisely

# A99-7

## FEATURES AND HIGHLIGHTS



### Technical parameters

▪ External (LxW):	1980x720mm
▪ Table height:	630-880mm
▪ Back plate up/down:	0-55°/0-18°
▪ Seat plate up/down:	0-35°/0-5°
▪ Leg board outwards:	0-90°

### Technical configuration


<input checked="" type="checkbox"/> Bed Frame	1 pc
<input checked="" type="checkbox"/> Foot Pedal	1 pc
<input checked="" type="checkbox"/> Hand controller	1 pc
<input checked="" type="checkbox"/> Anesthesia Screen	1 pc
<input checked="" type="checkbox"/> Hand Bar	2 pcs
<input checked="" type="checkbox"/> Leg Support with Cushion	2 pcs
<input checked="" type="checkbox"/> Arm Support with Cushion	2 pcs
<input checked="" type="checkbox"/> Stainless Steel Sewage Basin	1 pc
<input checked="" type="checkbox"/> Caster	1 pc

# Accessories

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



Pull Handle



Armrest




Knee Curutch



Pedal Remote Controller



Stainless Steel Basin



Anesthesia Frame



Height adjustment

630-880<sub>mm</sub>



Back-rest up

55°



Back-rest down

18°



Seat plate up/down

0-35°/0-5°



Reverse trendelenburg

0-90°



## QUALITY CONTROL


### Material

The base seat is plastic cover and height adjustment system is made of iron, environmental plastic can be sterilized in 100°C, tensile upto 30MPa.



### Functions

The main bed and auxiliary bed are connected with each other, auxiliary table is easy for rotating and dismantling.




### Handle controller

Handle controller which is easy for moving obstetric table.




### Leg support

Double using leg support, can be used as boatstretcher.




### Mattress

The mattress adopts high density memory sponge and is once formed, with no seam.



### Bed width

seamless foam mattress, with waterproof PU cover. Easy cleaning and disinfection.




### Control system

Good quality domestic motor, Various actions are controlled by remote controller.



### Sewage basin

Equipped with stainless steel sewage basin and other various accessories.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### South America Branch

#### Saikang Medical Technology, Inc

W3400 Corrientes, Argentina  
 TEL: +55 (11) 2758-0438  
 Rua Universal, Hollywood,  
 São Bernardo do Campo,  
 Brazil  
 TEL: +55 (11) 2758-0438

### East Africa Branch

#### Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,  
 upper hill, nairobi, Kenya  
 TEL: +254714084806

### South East Asia Branch

#### Saikang Medical Technology Ltd

Sala Thammasop Subdistrict  
 Thawi  
 Watthana District  
 Bangkok Thailand  
 TEL: +66 642741046

### UK Branch

#### Uk Healthward International Limited

DEPT 111,  
 196 High Road Wood Green,  
 London, N22 8HH  
 United Kingdom  
 TEL: +31 (0) 202111106

### Russia Branch

#### LLC 'SAIKANG'

454004, Russian Federation,  
 Chelyabinsk,  
 Universitetskaya  
 Naberezhnaya, Russia  
 TEL: +7 (919)300-00-06

### Mid East Branch

#### Saikang Medical Supplies (Mid East) LLC

Musaffah  
 P.O Box 53200, Abu Dhabi,  
 UAE  
 Tel: +971 585710652

**SKSurgitech™**  
 by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
 STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ export@sksurgitech.com

🌐 www.sksurgitech.com



Company video



General catalogue