

# A99-7


## MULTI-FUNCTIONAL ELECTRIC OBSTETRIC BED

Registration code: SKF-C



# A99-7

## FEATURES AND HIGHLIGHTS



### Technical parameters

▪ External Size (LxWxH)	1980x720x(630-880) mm
▪ Back Plate Up/Down	55°/18°
▪ Seat Plate Up/Down	35°/5°
▪ Leg Plate Swing Out	90°



### Technical configuration

<input checked="" type="checkbox"/> Bed Frame	1pc
<input checked="" type="checkbox"/> Foot Pedal	1pc
<input checked="" type="checkbox"/> Hand Controller	1pc
<input checked="" type="checkbox"/> Anesthesia Screen	1pc
<input checked="" type="checkbox"/> Hand Bar	2pcs
<input checked="" type="checkbox"/> Leg Support with Cushion	2pcs
<input checked="" type="checkbox"/> Arm Support with Cushion	2pcs
<input checked="" type="checkbox"/> Stainless Steel Sewage Basin	1pc
<input checked="" type="checkbox"/> Caster	1pc



**Material**

The base seat is plastic cover and height adjustment system is made of iron, environmental plastic can be sterilized in 100°C, tensile upto 30MPa.



**Function**

The main bed and auxiliary bed are connected with each other, auxiliary table is easy for rotating and dismantling.



**Control System**

Handle controller which is easy for moving obstetric table.



**Leg Support**

Double using leg support , can be used as boatstretcher.



**Mattress**

The mattress adopts high density memory sponge and is once formed,with no seam.



**Stainless Steel Sewage Basin**

Equipped with stainless steel sewage basin and other various accessories.



**Control System**

Good quality domestic motor,Various actions are controlled by remote controller.



**Pull Handle**



**Armrest**



**Pedal Remote Controller**






## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🕒 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue