



A99-8

GYNECOLOGICAL EXAM COUCH

Registration code: SKL-B






SKSurgitech[™]
by Saikang Medical

Care you precisely

A99-8

FEATURES AND HIGHLIGHTS




Technical parameters

External Size (LxWxH)	1400x740x610-1100 mm
Back Plate Up/Down	70°/20°
Seat board upward	40°
Main Voltage AC	220V±22V,50Hz±1Hz
Patient Weight Capacity	185kg

Technical configuration

<input checked="" type="checkbox"/> Pull Handle	2pcs
<input checked="" type="checkbox"/> Knee Curtch	2pcs
<input checked="" type="checkbox"/> Pedal Remote Controller	1pc
<input checked="" type="checkbox"/> Stainless Steel Basin	1pc
<input checked="" type="checkbox"/> Instrution Manual	1pc



Material

The base seat and height adjustment system are made of powder coated steel, high gloss, anti corrosion, easy to clean.



Function

Double using leg support , can be used as boatstretcher.



Pedal Remote Controller

Famous brand motor,Various actions are controlled by remote controller.



Pull Handle

The hand bar can be dismantled and movable.



Stainless Steel Basin

Equipped with stainless steel sewage basin and other various accessories.



Mattress

The mattress adopts high density memory sponge and is once formed,with no seam.




Wheel

Each wheel equipped with single brake which is easy for moving obstetric table.



Knee Curtch





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa





Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 ·  24/7
 export@sksurgitech.com
 www.sksurgitech.com



Company video



General catalogue