

SDL-D0301

PORTABLE PRESSURE
STEAM STERILIZER



SDL-D0301

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Working pressure	0.14-0.16Mpa
▪ Working temperature	126°C
▪ Temperature adjust range	105~126°C
▪ Timer range	0~99min
▪ Max. Safety pressure	0.165Mpa
▪ Heat average	≤±1°C
▪ Power	AC220V.50Hz/2KW

Characteristics

- Fully stainless steel structure.
- Automatically shut off with beep reminding after sterilization
- Hand wheel type of quick-open door structure.
- Electric heated.
- Digital display of working status, touch type key
- Easy to operate safe and reliable.

Description

This product is a rapid and effective disinfection and sterilization equipment using pressure cell and steam. It is suitable for hospitals, health centers, health stations, infirmary of factories and mines, scientific research units and so on. It can sterilize medical instruments, efficients, glassware, bath liquid, culture medium, etc. It can also be used as cooking utensils and high quality drinking water in plateau area.



The material
High quality stainless steel material



Digital display
LED display automation





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue