

SDL-D0313

HORIZONTAL PRESSURE STEAM STERILIZER



SDL-D0313

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Working pressure	0.22Mpa
▪ Working temperature	134°C
▪ Temperature adjust range	40-134°C
▪ Time for sterilization	0-60min
▪ Time for drying	0-60min
▪ Heat average	≤±2°C
▪ Power	AC380V.50Hz/9KW

Characteristics

- Automatically controlled the sterilizing process, easy to operate.
- The safety valve will be open automatically when the inner pressure over 0.24MPa and the steam would be exhausted the water tank.
- With drying function,suitable for the medical dressing drying
- With over-temperature,over-pressure auto-protect
- The door can not be open until the pressure in the chamber is reduced to 0.027MPa.
- And the unit can not be started on if the door doesn't close well.
- Automatically cut off the power if the water-lack happened and alarm at the meanwhile
- The chamber of the sterilizer is made of stainless steel

Description

This horizontal pressure sterilizer. Gravity displacement method is adopted to eliminate cold air in the inner chamber, which makes the content of cold air in the inner chamber less and the cause of destruction more manageable. At the same time, in the process of destroying the equipment, the control system automatically adjusts the inlet and outlet steam according to the internal chamber degree. It is mainly used in hospital supply room, operating room, pharmaceutical factory and laboratory of scientific research institute to sterilize cloth articles, surgical instruments and instruments used for medical and experimental purposes.



Drying function

With drying function, to achieve the drying effect of sterilized articles.



Overcurrent protection

The power breaker automatically cuts off the power when the current is overloaded.



Protection device

Can not start until the water level is high, low water level is short of automatic power off protection.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ export@sksurgitech.com

🌐 www.sksurgitech.com



Company video



General catalogue