

SDL-D0364B

TABLE TOP STEAM STERILIZER



SDL-D0364B

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Working pressure	0.22Mpa
▪ Working temperature	134°C
▪ Temperature adjust range	105-134°C
▪ Timer range	0~99min
▪ Heat average	≤±1°C
▪ Power	AC220V.50Hz/2.7KW

Characteristics

- For 4~6 minutes rapidly sterilizing.
- Safe protection of water lacking
- Digital display of working status, touch type key
- Exhaust the cool air automatically
- With 3 fixed cycles of sterilizing program and adjustable for user
- Door safety lock system.
- The whole process of water adding, temperature rising, sterilizing, drying steam discharging controlled automatically.
- With four stainless steel sterilizing plates.
- The chamber of the sterilizer is made of stainless steel.
- Steam-water inner circulation system: no steam discharge, and the environment for sterilizing will be clean and dry
- Automatically shut off with beep reminding after sterilization
- With drying function.

Description

The sterilizer is a safe, reliable and automatic control of small rapid desktop sterilization equipment, suitable for medical health, pharmaceutical, scientific research and other units. For surgery, dentistry, ophthalmic instruments. Glassware, pharmaceuticals, media and fabric dressings, and foods that are resistant to saturated steam provide rapid sterilization.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue