

SK-D001

DEFIBRILLATOR



SK-D001

FEATURES AND HIGHLIGHTS

Optional

Defibrillation pulse: Monophasic/Biphasic

Appearance

The whole defibrillator is small and portable, novel design, unique and beautiful shape, with good control.




Paddle

Composite paddles, reusable external adult paddles, pediatric paddles integrated. Text instructions at the handle for easy operation by medical staff. Discharge by pressing buttons on paddles together.



Energy

Energy step: 0, 3, 5, 7, 10, 20, 30, 50, 100, 200, 300, 360 Joules.
Charging time: at 360J less than 9s; at 100J less than 2s.



Battery

Power supply: mains(AC100~240V) and/or battery.
Battery: sealed rechargeable lead-acid battery(optional lithium battery).
Battery charge 3.5h to full, full charge to 35 shocks(360J).



Back

1. The AC switches, 'I' state indicates power on, 'O' state indicates power off.
2. Fuse specs: T3AL250V, $\Phi 5 \times 20$ mm.
3. Potential equalization conductor.
4. AC input socket.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue