



# SK-EH003

OXYGEN-BOTTLE



SKSurgitech™  
by Saikang Medical

Care you precisely

# SK-EH003 Oxygen-bottle

## Technical parameters

▪ Empty bottle weight:	6 kg
▪ Water volume:	4 l
▪ Thickness:	2.8 mm
▪ Content:	Oxygen
▪ Hydraulic test pressure:	22.5 MPa
▪ Nominal working pressure:	15.0 MPa
▪ Air tightness test pressure:	15.0 Mpa

## Technical configuration

<input checked="" type="checkbox"/> Inflatable oxygen bridge	1 pc.
<input checked="" type="checkbox"/> Humidifying bottle	1 pc.
<input checked="" type="checkbox"/> Oxygen tube	1 pc.
<input checked="" type="checkbox"/> Wrench	1 pc.
<input checked="" type="checkbox"/> Case	1 pc.



### Best quality materials

The bottle material is 37 Mn steel.  
Standard: GT/T5099.3-2017  
(Steel Seamless Bottle)



### Long life time

Each bottle has a unique steel engraved code. The life time of bottle (accessories not included) is 30 years.



**SKSurgitech**<sup>TM</sup>  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎ +86 18021231901 · 🕒 24/7  
✉ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue