



SK-EH303

AV-2000B3 VENTILATOR



SKSurgitech[™]
by Saikang Medical

Care you precisely

SK-EH303

FEATURES AND HIGHLIGHTS

The characteristics of

- Display mode 10.4-inch High-definition TFT color LCD screen display
- Ventilation mode VCV,PCV,SIMV.PSV.PSV+SIMV,CPAP
- Ventilation function PEEP, SIGH, standby, IP, IRV
- Ventilation parameters
Tidal volume 20~1800ml
Rate 2~120bpm
SIMV rate 2~20bpm
I:E 4:1~1:8
Pressure trigger sensitivity -10~10cmH2O
Flow trigger sensitivity 1~20L/min
PEEP 0~25cmH2O
Pressure range 3~50cmH2O
Pressure limitation 5~60cmH2O
Pressure slope 1~10gear
Flow rate setting 1~10gear
Expiratory trigger sensitivity 10%~90%
Time control 1~10s
SIGH 1.5 times the inspiratory time
(60~120adjustable)
Inspiratory platform 0~50%
Oxygen concentration 21%~99%



Technical parameters

- Parameters for ventilation monitoring
Expired tidal volume, Inhale tidal volume, Minute Ventilation volume, IPPV rate, SIMV rate, total respiratory rate, autonomous breathing rate/Epeak pressure of airway average pressure, PEEP, inspiratory trigger pressure, inspiratory platform, oxygen concentration, battery capacity, pressure-time waveform, flow rate-time waveform, flow-volume loop, pressure-volume loop
- Security Alarm System
Oxygen concentration alarm upper limit setting range 21%~99%
low limit setting range 18%~80%
Airway pressure alarm upper limit setting range 5~60cmH2O
low limit setting range 0~50cmH2O
minute ventilation volume alarm upper limit setting range 1~20L/min
low limit setting range 0~19L/min
Tidal volume upper limit setting range 10~1800ml
low limit setting range 0~1700ml
Respiratory rate upper limit setting range 10~99bpm
low limit setting range 0~50bpm
Sustained high-pressure alarm it will give alarm when stress have consistently been higher than 30cmH2O
Suffocation alarm Power alarm
Low pressure of Oxygen alarm Intubation off alarm
Low pressure of Air alarm battery voltage alarm
- Power AC 220V 50Hz
UPS, Storage battery
- Mechanical Arms
- Air compressor With "standby" function



CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎ +86 18021231901 · 🕒 24/7
✉ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue