

# SK-EM032

INTEGRATED  
HEALTHCARE  
MONITOR



# SK-EM032

## FEATURES AND HIGHLIGHTS

### Introduction

Integrated healthcare monitor is newly designed in accordance with clinical and hospital examination

### Descriptions

- 8.4" / 10.4" / 12.1" high resolution industry color TFT LCD screen.
- Basic 5 standard parameters: ECG, NIBP, SpO2, RESP, 2-Temp.
- 3/5 lead, max 7 ECG wave forms display.
- Humanized operation interface; Multi-language available.
- Audible & visual alarm and adjustable alarm limits.
- Mini USB/OTG port, SD card slot, SIM (3G)port, USB for transferring patient's data to PC.
- 120 hours full discloser and 1200 groups of data storage capacity, 500 groups of NIBP stored.
- Rechargeable high-energy li-ion battery.
- Arrhythmia analysis and S-T segment analysis



### Standard accessories



5-Lead Integrated Button Connector ECG Cable



Reusable Finger Clip SpO2 Senso (ADU)



NIBP Tube (for Reusable Cuff)



Skin Temperature Probe (10KΩ)



Reusable Adult Cuff



Power cord



## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi,  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai,  
United Arab Emirates.  
Tel.: +20 128 8962497

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🌐 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue