

SK-EM032

INTEGRATED
HEALTHCARE
MONITOR



SK-EM032

FEATURES AND HIGHLIGHTS

Introduction

Integrated healthcare monitor is newly designed in accordance with clinical and hospital examination

Descriptions

- 8.4" / 10.4" / 12.1" high resolution industry color TFT LCD screen.
- Basic 5 standard parameters: ECG, NIBP, SpO2, RESP, 2-Temp.
- 3/5 lead, max 7 ECG wave forms display.
- Humanized operation interface; Multi-language available.
- Audible & visual alarm and adjustable alarm limits.
- Mini USB/OTG port, SD card slot, SIM (3G)port, USB for transferring patient's data to PC.
- 120 hours full discloser and 1200 groups of data storage capacity, 500 groups of NIBP stored.
- Rechargeable high-energy li-ion battery.
- Arrhythmia analysis and S-T segment analysis



Standard accessories



5-Lead Integrated Button Connector ECG Cable



Reusable Finger Clip SpO2 Senso (ADU)



NIBP Tube (for Reusable Cuff)



Skin Temperature Probe (10KΩ)



Reusable Adult Cuff



Power cord



CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🕒 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue