

# SK-EM201

## SYRINGE PUMP



# SK-EM201

## FEATURES AND HIGHLIGHTS

**Accessories:** host, remote control, power cord, user manual, warranty card, product certificate.

**Application:** The main application software is developed on the Linux platform, which is more secure and reliable. The function and technology are completely independently researched and developed, and can fully meet customers' future upgrade needs.

**Other Functions:** Automatic syringe specification recognition function: can be used for 5ml, 10ml, 20ml, 30ml, 50 (60) ml syringes. The manually adjustable precision function is more convenient for adjusting and ensuring precision requirements



### Technical parameters

- Product weight gross weight 3kg net weight 2.3kg

---

- Injection modes speed mode, time mode, weight mode (2 types).

---

- Injection volume range 0.1~1999.9ml (in increments of 0.1ml)

---

- Cumulative volume range 9999.9ml (in increments of 0.1ml)

---

- Accuracy of injection speed  
within  $\pm 2\%$  (not including the error of the syringe itself)

---

- Injection blockage  
Three levels of blockage alarm pressure (high, medium, low) can be selected according to clinical needs

---

- KVO speed (Keep vessels open)  
0.1~5ml/h (adjustable in increments of 0.1 ml/h)

---

- Remote Control Function Distance range 0m--5m

---

- Working environment Temperature: +5°C~+40°C

---

- Relative humidity 20%~ 90%

---

- Storage environment Temperature -30°C~+55°C

---

- Atmospheric pressure 860hPa~1060hPa

---

- Safety level  
Class I CF equipment, built-in power supply, continuous operation with intermittent loading

---

#### Injection speed range: 50ml syringe:

0.1~1800ml/h (0.1~999.9ml/h increment is 0.1ml/h  
1000~1800ml/h increment is 1ml/h)  
30ml syringe: 0.1~900ml/h  
20ml syringe: 0.1~600ml/h  
10ml syringe: 0.1~300ml/h  
5ml syringe: 0.1~150ml/h

#### Fast forward speed:

50ml syringe: 1800ml/h (Bolus speed 1200ml/h)  
30ml syringe: 900ml/h (Bolus speed 600ml/h)  
20ml syringe: 600ml/h (Bolus speed 400ml/h)  
10ml syringe: 300ml/h (Bolus speed 200ml/h)  
5 ml syringe: 150ml/h (Bolus speed 100ml/h)

#### Alarm function:

alarm for injection completion, syringe empty, injection near, injection blockage, syringe falling off, syringe installation error, battery low (undervoltage alarm), AC power disconnection, stop timeout prompt.

**Design**

Humanized integrated design, exquisite appearance, precise injection, and powerful functions.



**Shell Material**

New ABS engineering plastic injection molding.



**Screen Design**

4.3-inch LCD screen, support Chinese and English language. Equipped with backlight, it can be used in various lighting conditions. The operating interface is presented at a 45-degree angle, which enhances the clarity of the visual effect.



**Internal Structure**

The aluminum-magnesium alloy bracket is stable and reliable, has good shock resistance and is easy to move.



**Handle**

Fixed handle design, comfortable, stable and easy to move.



**Syringe Fixture**

Using an open-and-close design, it is easy to clamp the syringe, which can accommodate a 60mm syringe and is more in line with clinical operation.



**Operation Interface**

Interface display content: time, battery capacity, injection status, mode, injection speed, injection volume, cumulative volume, injection time, syringe specification, alarm volume, block pressure, precision, weight, medication amount, dose, and liquid volume.



**Remote Control Function**

It can remotely control devices, directly input numbers to adjust parameters such as speed, infusion time, and medicine volume, which can save time and effort for medical staff.



**Clock Display Function**

It is convenient to memorize and check time.



**Alarm Function**

The triple alarm system of sound, light and screen, three levels of adjustable volume.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi,  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai,  
United Arab Emirates.  
Tel.: +20 128 8962497

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🌐 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue