

SK-EX101

ELECTRIC PORTABLE SUCTION MACHINE



SK-EX101

FEATURES AND HIGHLIGHTS

Specifications

▪ Model	DFX-23B+I
▪ Power supply	220V±10%50Hz+2%
▪ Air pumping efficiency	≥28L/min
▪ Adjustment range of negative pressure	0.013~0.09MPa
▪ Noise level	≤55dB
▪ Input power	150VA
▪ Liquid vessel	1000mlx1

Features

- AC and DC power supplies can be connected to car-carried power supply
- Self-storage supply up to 2 hours
- AC and DC automatic switching does not affect normal use imported membrane pump with longer life
- With full charge function, when the device is full charged the charging will be automatically stopped.
- Small size, especially for ambulance and field application
- With overflow protection unit to prevent liquid from getting into pump body
- Adopts high quality cold plate shell
- Surface electrostatic spraying
- Portable





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue