



SK-EX102

MOBILE ELECTRIC
SUCTION
DEVICE



SKSurgitech[™]
by Saikang Medical

Care you precisely

SK-EX102

FEATURES AND HIGHLIGHTS



TECHNICAL PARAMETERS

- Power voltage: AC220V±10% 50Hz
- Negative pressure: $\geq 0.09\text{MPa}$ (680mmHg)
- Suction rate: $\geq 20\text{L}/\text{min}$ (35L pump suitable)
- Range of negative pressure: 0.013MPa~0.09Mpa(680mmHg)
- Suction bottle: 2500ml x 2
- Noise: $\leq 55\text{dB}$
- Input power: 400VA
- Pump structure: Diaphragm type
- Working type: intermittent loading, continuous operation
- All - plastic ABS enclosure that is formed in one piece is adopted
- This device is equipped with a imported membrane pump available with long service life
- This device has low noise and is maintenance - free
- This device is equipped with an overflow protection unit to prevent liquid from getting into pump body
- It's equipped with the hand-operated switch and footstep - operated switch

PUMP
Imported membrane pump for long service life

A detailed cutaway view of a blue and white membrane pump mechanism, showing internal components like the diaphragm and motor.

DICL
Large dial data display for easy reading

A circular pressure gauge with a black dial and white markings. The needle points to approximately 100 mmHg. The dial has 'MPa' and 'mmHg' scales. The text 'Intermittent operation' and 'PRESSURE' are visible on the gauge's frame.

BOTTLE
5L bottle, suitable for big flow usage requirements

Two clear plastic bottles with black caps and white tubing attached, sitting on a surface.

SWITCH
have manual switch and foot switch, easy to operate

A grey foot pedal with a black cable and a manual switch on top.

FILTER
Filters improve hygiene standards

A collection of small, white, cylindrical filters with blue caps, arranged in a row.

OUTWARD
One-time molding of outer ABS case, seamless and durable

A full view of the white SK-EX102 device with two bottles.

PROTECTIVE DEVICE
Overflow protection device prevents liquid from entering the pump

A cylindrical metal component with a threaded end, used for overflow protection.



CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ export@sksurgitech.com

🌐 www.sksurgitech.com



Company video



General catalogue