

SK-EX103

ELECTRIC SUCTION MACHINE



SK-EX103

FEATURES AND HIGHLIGHTS

Features

- 23C.V: all plastic enclosure which is formed in one piece 23C.Jll adopts high quality cold plate shell
- Imported membrane pump with longer life
- Low noise and maintenance-free
- With overflow protection unit to prevent liquid from getting into pump body
- Be equipped with hand switch and foot switch

Specifications

▪ Model	DFX-23C-111
▪ Power supply	220V+10% 50Hz+2%
▪ Air pumping efficiency	≥40L/min 80L/min
▪ Adjustment range of negative pressure	0.013~0.09MPa
▪ Noise level	≤65dB
▪ Input power	4DOVA
▪ Liquid vessel	2500mlx2



Introduce

- (1) Imported diaphragm pump
- (2) Bottle is easy for cleaning & disinfection
- (3) Oil-free pump is maintenance-free



Introduce

- (4) Low noise
- (5) With an overflow protection
- (7) Working type: intermittent loading, continuous operation





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa






Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 ·  24/7
 export@sksurgitech.com
 www.sksurgitech.com



Company video



General catalogue