

SK-EX114

ELECTRIC SPUTUM
SUCTION MACHINE



SK-EX114

FEATURES AND HIGHLIGHTS

Specifications

▪ Model	DFX-23A-1
▪ Power supply	220V±10% 50Hz±2%
▪ Air pumping efficiency	≥15L/min
▪ Adjustment range of negative pressure	0.013~0.08MPa
▪ Noise level	≤65dB
▪ Input power	150VA
▪ Liquid vessel	1000ml

Features

- All plastic ABS shell that is formed in one piece(23.J23A.11)
- Macromolecule suction bottle is adopted. easy to disinfect and clean
- Oil-free pump and maintenance-free
- Small size and low noise
- With overflow protection unit to prevent liquid from getting into pump body





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🕒 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue