

# SK-EX128

ELECTRIC SPUTUM  
SUCTION MACHINE



# SK-EX128

## FEATURES AND HIGHLIGHTS

### Specifications

▪ Model	DFX-23A-1
▪ Power supply	220V±10% 50Hz±2%
▪ Air pumping efficiency	≥15L/min
▪ Adjustment range of negative pressure	0.013~0.08MPa
▪ Noise level	≤65dB
▪ Input power	150VA
▪ Liquid vessel	1000ml

### Features

- All plastic ABS shell that is formed in one piece(23.J23A.11)
- Macromolecule suction bottle is adopted. easy to disinfect and clean
- Oil-free pump and maintenance-free
- Small size and low noise
- With overflow protection unit to prevent liquid from getting into pump body





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi,  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai,  
United Arab Emirates.  
Tel.: +20 128 8962497

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🌐 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue