



SK-L002

35W HALOGEN MEDICAL
EXAMINATION LIGHT



SKSurgitech™
by Saikang Medical

Care you precisely

SK-L002

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Voltage	AC90~240V 50Hz-60Hz
▪ Power of bulb	halogen 35W
▪ Color temperature	4000K
▪ Life of bulb	2000h
▪ Illumination(200mm)	≥40000Lux
▪ N.W.	1.5KG
▪ Package dimension (with clamp)	600x120x260mm
▪ Package dimension (with floor stand)	845x215x155mm
▪ Packaging material	Carton

Description

It adopt 12V35W halogen bulb as light source. The lamp head can be rotated from multiple angles, the lamp pole can be bent freely, and the lighting from any angle can be realized by arbitrary positioning. There are table type, guide rail type, wall type, mobile floor type optional.



Application

Clinic, Gynecological examination, Pet medical, Dental examination, Beauty, Dental, Aesthetic Medicine, Product testing, etc.

Handle design

Avoid contact with high temperature lighthouse, humanized.



Hose bent

Multi-angle lighting.



Universal Joint

Multi-angle lighting.



Multiple fixed methods

With the floor frame, two kinds of lifting and fixed options.




Different installation methods

With wall fixture, guide rail fixture, desktop fixture. Three different installation methods.



Halogen light source

Small size, economical and practical.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue