



# SK-L009

9W LED MEDICAL  
EXAMINATION LIGHT




# SK-L009

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ Voltage	AC90~240V 50Hz-60Hz
▪ Power of bulb	LED 3Wx3
▪ Color temperature	4000-5000K
▪ Life of bulb	≥50000h
▪ Illumination	8000Lux(400mm)
▪ N.W.	1.5KG
▪ Package dimension (with clamp)	600x120x260mm
▪ Package dimension (with floor stand)	845x215x155mm
▪ Packaging material	Carton



### Description


It adopt 3pcs 3W international brand LED bulb as light source,with high brightness,no temperature rise.  
 Brightness adjustable,the flexible hose bends freely for arbitrary positioning and multi-angle lighting.  
 There are table type, guide rail type, wall type, mobile floor type optional.

### Application

Clinic, Gynecological examination, Pet medical, Beauty, Dental, Product testing, Jewelry identification, etc.

#### Gooseneck design

Easier to locate and hover without sagging.



#### Flood lighting

More lighting area means less shadow.




#### Brightness adjustable



#### Multiple fixed methods

With the floor frame, two kinds of lifting and fixed options.



#### Different installation methods

With wall fixture, guide rail fixture, desktop fixture. Three different installation methods.



#### LED light

Higher brightness, lower consumption, longer service time.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi,  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai,  
United Arab Emirates.  
Tel.: +20 128 8962497

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🕒 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue