



SK-L009

MEDICAL EXAMINATION LIGHT



SKSurgitech[™]
by Saikang Medical

Care you precisely

SK-L009

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Voltage	AC90~240V 50Hz-60Hz
▪ Power of bulb	LED 3Wx3
▪ Color temperature	4000-5000K
▪ Life of bulb	≥50000h
▪ Illumination	8000Lux(400mm)
▪ N.W.	1.5KG
▪ Package dimension (with clamp)	600x120x260mm
▪ Package dimension (with floor stand)	845x215x155mm
▪ Packaging material	Carton



Description

It adopt 3pcs 3W international brand LED bulb as light source,with high brightness,no temperature rise.
 Brightness adjustable,the flexible hose bends freely for arbitrary positioning and multi-angle lighting.
 There are table type, guide rail type, wall type, mobile floor type optional.

Application

Clinic, Gynecological examination, Pet medical, Beauty, Dental, Product testing, Jewelry identification, etc.


Lamp head adjustment

The lamp head angle is easy to adjust and the angle of adjustment is wide.



Bulb

It uses imported bulbs, the average service life is up to 50000 hours.



Base

Optional base with swivel casters, easy to move.



Height adjustment

Height is adjustable, can match the height of the patient easily.



Different installation methods

Quick changing of clip type, wall type and floor type.



Color temperature&Illumination

Using LED light source, the illuminance can reach 50000 Lux, the color temperature is around 5500K, light spot and brightness are adjustable.





 **CONSTANT
INNOVATION**



Design



Research



Manufacturing



Maintenance



Service

 **WORLDWIDE
BRANCHES**

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa




Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 ·  24/7
 export@sksurgitech.com
 www.sksurgitech.com



Company video



General catalogue