



SK-L027

21W LED MOBILE MEDICAL
EXAMINATION LIGHT



SKSurgitech™
by Saikang Medical

Care you precisely

SK-L027


FEATURES AND HIGHLIGHTS

Technical parameters

▪ Voltage	AC90~240V 50Hz-60Hz
▪ Power of bulb	LED 3Wx7
▪ Color temperature	4000-5000K
▪ Life of bulb	≥50000h
▪ Illumination	18000Lux
▪ N.W.	6.5KG
▪ Package dimension	1220x260x240mm
▪ Packaging material	Carton

Description

It adopt 7pcs 3W international brand LED bulb as light source,with high brightness,long life span.
There is a big handle on lamp head, switch in front, make the lamp easy to operate,flexible hose could achieve multi-angle lighting.




Application

Gynecological examination, General diagnosis and treatment, Dental examination, Cosmetology, Pet medical, etc.

Gooseneck design

Easier to locate and hover without sagging.



Flood lighting

More lighting area means less shadow.



Wide voltage power supply

95V-240V wide voltage power supply compatible with all countries.



Special support



Different installation methods

With wall fixture, guide rail fixture, desktop fixture. Three different installation methods.



LED light

Higher brightness, lower consumption, longer service time.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ export@sksurgitech.com

🌐 www.sksurgitech.com



Company video



General catalogue