



SK-L309

WALL-MOUNTED
SURGICAL LIGHT



SKSurgitech[™]
by Saikang Medical

Care you precisely

SK-L309

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Voltage	AC90~240V 50Hz~60Hz
▪ Bulb power	LED 3Wx36
▪ Color temperature	4000~5000K
▪ Bulb life	≥ 50000h
▪ Working distance: 1000mm	
▪ Center illumination	≥65000Lux
▪ CRI	85≤Ra≤100
▪ Minimum spot diameter	≤125mm
▪ Lighting depth	≥200mm
▪ Temperature rise	≤2°C
▪ Lamp holder diameter	500mm
▪ Maximum rotation radiu	1800mm
▪ Installation height (recommended)	2500mm
▪ Wall requirements	concrete structure, brick structure, 240mm brick body, steel plate structure
▪ Net weight	28KG
▪ Packing size	1060*600*390mmmm
▪ Packing material	wooden case



Description

This product adopts 36 pieces 3W LED bulbs as lght source. High and adjustable brightness,small lighting spot.Lighting depth is up to 200mm.

Base castor configuration, easy to move;Handle control,lamp holder can rotate freely get multi-angle lighting

Application

Outpatient surgery, Oral implants, Plastic surgery, Pet medical treatment

Lamp panel

Dimeter of the lamp panel is 500mm



Floodlighting

The maximum exposure area means the minimum shadow



Brightness

Brightness is adjustable



LED technology

Higher brightness,lower consuivption,longer service time



Optional

Wall mounted type and Ceiling type





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue