

# SK-L403

AC/DC HEADLIGH



# SK-L403

## FEATURES AND HIGHLIGHTS


### Technical parameters

▪ Voltage	AC 90-240V 50-60Hz
▪ Bulb Power	LED 3W
▪ Color temperature	4000~5000K
▪ Bulb life	≥50000h
▪ Illumination	≥12000Lux
▪ DC charging time	≤5h
▪ Power supply time	≥4h
▪ Load weight on head	200g
▪ Net weight	450g
▪ Package size	230x190x130mm
▪ Package material	paper box

### Description

This product adopt 3W high-power international brand LED as light source, By optical design, has uniform light distribution, high brightness, light spot is clear and size adjustable.

This headlight could work while charging.



### Application

ENT, Dental, General surgery, Plastic surgery, Pet medical.

#### Adjustable light spot

The spot size can be adjusted freely according to the demand.



#### Lithium battery

AC and DC dual-purpose, can work while charging.




#### Adjustable lighthouse angle

The lighthouse angle can be adjusted freely according to the demand.



#### Magnifier installation joint





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi,  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai,  
United Arab Emirates.  
Tel.: +20 128 8962497

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🌐 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue