

SK-L406

5W HEADLIGHT
WITH HIGH-BRIGHTNESS



SK-L406

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Voltage	AC 90-240V 50-60Hz
▪ Bulb Power	LED 5W
▪ Color temperature	4000~5000K
▪ Bulb life	≥50000h
▪ Illumination	≥18000Lux
▪ DC charging time	≤5h
▪ Power supply time	≥4h
▪ Load weight on head	≤230g
▪ Net weight	470g
▪ Package size	230x230x140mm
▪ Package material	paper box

Description

It adopt 5W LED as the light source,with high brightness and long life span.

The light distribution is uniform,light spot clear and bright,brightness adjustable.



Application

ENT, Dental, Micro-plastic, Pet medical, Clinic examination.

Features

Physical heat dissipation.
No low frequency interference.
Long life.



Ergonomic design

The headband has multiple adjusting points and a detachable forehead pad.



Adjustable lighthouse angle

The lighthouse angle can be adjusted freely according to the demand.



Magnifying glass

It can be used with magnifying glass.



High power headlamp

5W high power LED output.



High efficiency lithium battery

It can with continuously for up to 4 hours.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue