

# SKH007-1

## INSTRUMENT TROLLEY

Registration code: SKM-B



# SKH007-1

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ Size L (LxWxH)	1000×500×900 mm
▪ Size M (LxWxH)	880×430×860 mm
▪ Size S (LxWxH)	760×360×820 mm

### Technical configuration

<input checked="" type="checkbox"/> 3" Castor	4pcs
---	------



#### Application

Stainless steel instrument trolley, Multifunctional nursing trolley, Medical trolley.



#### Material

Φ25mm stainless steel round tube, 0.8mm thick stainless steel board. The whole frame used high quality S.S. with the mirror surface polishing.



#### Craftsmanship

Imported intelligent robot is used for welding forming, high precision, no blind spots, smooth surface, no burrs, beautiful and strong, high bearing capacity.



#### Triple

Triple stainless steel apply mechanically operative instrument table, without siderail, the smallest one with two-layer countertop design.



#### Castors

3 inch covered castors (diagonal brake), flexible and light to drive, smooth steering, wear-resistant, two with brake, can be easily braked in any state.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi,  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai,  
United Arab Emirates.  
Tel.: +20 128 8962497

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🕒 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue