



# SKH038-6

## STAINLESS STEEL TROLLEY

Registration code: SKM-B



**SKSurgitech**<sup>™</sup>  
by Saikang Medical

Care you precisely

# SKH038-6

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ Overall dimension	600×400×800-1100mm
▪ Tray size	600×400×20mm
▪ Safe load-bearing	≥5KG
▪ Product weight	9KG



#### Design

Sophisticated design, easy to install, stable and durable.



#### Material

High-quality stainless steel, better anti-corrosion and not easy to rust.



#### Tray

Stainless steel tray thickness:0.8mm, size: 600\*400\*20mm, raised edges around the tray to prevent the falling of items.




#### Column

Stainless steel column diameter: Ø60mm, thickness: 1.2mm. Height can be adjusted hydraulically and manually by foot pedal, adjustable height range: 800-1100mm. Two optional lifting modes, easy and convenient.



#### Base

Stainless steel base diameter: Ø25mm, thickness: 1.2mm, U-shaped design for excellent stability.



#### Foot pedal

The wide-surface and antislip design ensures a comfortable and secure footing.



#### Castor

4 pcs anti-static castors with brake function.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi,  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai,  
United Arab Emirates.  
Tel.: +20 128 8962497

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🌐 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue