

SKH048

OXYGEN BOTTLE TROLLEY

Registration code: SKM-B



SKH048

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 1435×360×1040 mm


Technical configuration

- ✓ 6" TPP Rubber Rastor 2pcs
- ✓ 3" Covered Castors 2pcs



Design

U-shaped groove design, the radian meets the requirements of the oxygen bottle for safe transportation.



Steels | Welding

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth wedding.



Dimension

Size: 1435×360×1040mm, Can load 40L oxygen cylinders, the body is light and easy to move.



Galvanized Chain

4mm galvanized chain, operational safety.



Handle

Pushing handle, easy to move the trolley.



Thickening Plate

Cold rolled steel plate, effective compression, bearing weight, not easy to damage, long service life.



6" TPP Rubber Rastor

6" TPP rubber castor, Strong, durable, labor-saving, safty in use.



3" Covered Castors

3" Covered Castors, each with brakes, easy to install, stable and smooth.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue