



SKH070

WORK TABLE

Registration code: SKM-B



SKSurgitech[™]
by Saikang Medical

Care you precisely

SKH070

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 1800×700×800 mm



Material

The whole frame used high quality S.S. with the mirror surface polishing, easy to clean and wipe, not easy to rust.



Craftsmanship

Imported intelligent robot is used for welding forming, high precision, no blind spots, smooth surface, no burrs, beautiful and strong, high bearing capacity.




Design

Double-layer work table with stainless steel plate for the work surface and rectangular stainless steel tube for the second layer.



Table foot

Work table feet with high quality rubber sleeve, avoid direct contact between the table foot and ground to damage the ground.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue