



SKN001

STAINLESS STEEL TANK



SKN001

FEATURES AND HIGHLIGHTS

Containers used for storage of sterilization accessories or surgical instruments to be sterilized by autoclaving.



5 kinds of specifications are available to meet different needs.



Material

Strictly selected 304 stainless steel, thickened material, anti-rust, resistant to falling and not easy to break.



Design

Clamshell design to prevent dust. With handle, easy to carry and prevent from falling.

The inner wall is polished with fine sanding, leaving no dirt, making cleaning more convenient.



Hole

The edge is added with a ventilation hole switch lock design, and the ventilation holes can be adjusted freely.



Base

The base is designed with small feet to prevent slippage and scratches on the bottom, and has high wear resistance.



CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue