



# SKR-IB00

## SIMPLE LAPTOP CART

Registration code: SKM-A



**SKSurgitech**<sup>™</sup>  
by Saikang Medical

Care you precisely

# SKR-IB00

## FEATURES AND HIGHLIGHTS

### Technical parameters

- Size (LxWxH) 533×455×750-1150 mm

### Technical configuration

<input checked="" type="checkbox"/> Big Basket	1pc
<input checked="" type="checkbox"/> Tray	1pc
<input checked="" type="checkbox"/> Castor	4pcs



#### Column

Aluminum alloy column, the top equipment mesa.



#### Table Top

ABS material top, easy to clean and durable to scratch. Cover is free to open, laptop can be put in easily.



#### Compartement

Compartement can storage any laptop size under 17", is underneath top cover and lockable with key.



#### Holder for Mouse,

Double sided telescopic holder for mouse, left or right hand operating is both available.



#### Base

Base is ABS covered high quality metal structure, featured with both stability and easy-maintenance.



#### Castor

Casters is using medical silent material, double sided casters, stable and noiseless. Front 2 casters with brake, lockable at 360 degree, rear 2 casters free and swivel to any direction.



#### Big Basket

With basket on the back, easy to place goods, Environmentally plastic is non-toxic and easy to clean.



#### Option

Option configuration of high performance li-battery system, with LED Keyboard light, constantly providing stable power source, safe and convinient.






## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SKSurgitech™

by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🕒 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue