



# SKR-J10

## PATIENT MONITOR CART

Registration code: SKM-A



**SKSurgitech**<sup>™</sup>  
by Saikang Medical

Care you precisely

# SKR-J10

## FEATURES AND HIGHLIGHTS

### Technical parameters

- Size (LxWxH) 505×430×950 mm

### Technical configuration

- Big Basket 1pc
- Tray 1pc
- Castor 4pcs



#### Column

Aluminum alloy column, the top equipment mesa.



#### Computer Holder

Concise and compact for 13" patient monitor, Proposal Table load  $\leq 15$ kg.



#### Big Basket

With basket on the back, easy to place goods, Environmentally plastic is non-toxic and easy to clean.



#### Tray

Built-in ABS tray, ABS tray size: 388mm (L) × 238mm (W) × 50mm (H), the height of the ABS tray above the ground is 579.8mm.



#### Base

Base is ABS covered high quality metal structure, featured with both stability and easy-maintenance.



#### Castor

4pcs of 75mm double side casters, each with brake.






## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SKSurgitech™

by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🕒 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue