

# SKR-J21

## ENDOSCOPIC SYSTEM EQUIPMENT CART

Registration code: SKM-A



# SKR-J21

## FEATURES AND HIGHLIGHTS

### Technical parameters

- Size (LxWxH) 665×630×1650 mm

### Technical configuration

- Castor 4pcs



#### Column

Metal structure material comes from Baosteel from Fortune Global 500.




#### Holder of Screen

Above the table is equipped with a movable universal arm, which can hang the display, can rotate 360 degrees and adjust the height.



#### Working Table

Countertop height: 1125mm, dimensions: 527 (L) mm × 420mm (W), with a drawer under the countertop.



#### Movable Shelf

The height and number of each floor can be adjusted freely. The standard product is 2 adjustable movable shelves, which can be optionally added. The recommended load of the movable shelves is  $\leq 20$ kg.



#### Bumper Cover

With bumper, protect the cart during transportation.



#### Castor

4pcs of 100mm double side casters, each with brake.






## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SKSurgitech™

by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🕒 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue