



SKR-MT606

MEDICINE TROLLEY

Registration code: SKM-B



SKSurgitech™
by Saikang Medical

Care you precisely

SKR-MT606

FEATURES AND HIGHLIGHTS

Technical parameters

- External size (LxWxH) 1100x470x1000 mm
- Castors 125 mm



Material

The whole trolley is made of aluminum, plastic and steel, which is light and flexible, novel in design, unique and beautiful in shape, with good controllability.



Table top

Concave ABS tabletop, avoid possible falling down during movement. With push handles on two sides. The guardrail is hung with a detachable partition, which can be adjusted to place bottles of disinfectant cup, grinding wheel cup, etc.




Left Side

Detachable large white hanging box on the left: small instruments or gloves can be placed.



Hooks

Hanging small and light items; Edge box and net basket: place used needle cylinder and other articles.



Garbage can

It is convenient for medical staff to use, avoid contact infection.



Drawers

1-3 layer shallow drawers: 80mm high, 492x342x62 for internal separation;
4-6 layer deep drawers: 160mm high, 492x342x137 for internal separation, and oral medicine can be placed separately.



Castor

5 inch directional silent castors, flexible and light to drive, smooth steering, wear-resistant, two with brake, can be easily braked in any state.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue