




# SKR-R03

## LAPTOP CART

Registration code: SKM-A



SKSurgitech™  
by Saikang Medical



# SKR-R03

## FEATURES AND HIGHLIGHTS

### Product Introduction

The mobile doctor workstation based on wireless network integrates highly integrated computer, wireless communication, intelligent power management, mobile management technology, and is suitable for a variety of clinical work environments; whether it is outpatient, emergency, inpatient ward, operating room or ICU, it can help clinical Medical workers realize intelligent mobile ward rounds.

### Technical Parameters

▪ Dimensions L×W×H	530×520×850mm
▪ Work table L×W	445×410mm
▪ Sub-table size L×W	210×200mm
▪ Keyboard base size L×W	440×250mm
▪ Machine weight	21kg
▪ Lifting size	table lifting size 30mm
▪ Caster design	4 inch medical grade silent wheel



## Material

The table top and column are made of aviation-grade aluminum alloy, the base and handle are made of one-time molding all-cast aluminum material, the whole car adopts high-voltage electrostatic spraying process, and imported environmentally friendly antibacterial plastic powder is used to ensure that the product meets the medical and health requirements, and the surface is matte and non-dazzling.




## Table Top

Aviation grade aluminum alloy table, easy to clean and anti-corrosion, with partition design on both sides to prevent items from falling.




## Keyboard Subsurface

The pull-out sub-table can place the keyboard and documents, and provides a hidden extended mouse panel.



## Lower Column

The column is equipped with a professionally adjusted air rod, which is easy to rise and fall, and can be quickly adjusted to a suitable height, adapting to different sitting and standing postures, with a stroke of up to 20cm. And optional electric lift.



## Chassis

All-metal stretched X-shaped chassis, one-time cast aluminum molding, to ensure the center of gravity and level of the vehicle.



## Casters

4-inch dual-wheel medical grade silent wheels, all with brakes, after tens of thousands of kilometers of long-term durability tests, flexible and durable.





## CONSTANT INNOVATION

Design

Research

Manufacturing

Maintenance

Service



## WORLDWIDE BRANCHES

### South America Branch

#### Saikang Medical Technology, Inc

W3400 Corrientes, Argentina  
TEL: +55 (11) 2758-0438  
Rua Universal, Hollywood,  
São Bernardo do Campo,  
Brazil  
TEL: +55 (11) 2758-0438

### East Africa Branch

#### Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,  
upper hill, nairobi, Kenya  
TEL: +254714084806

### South East Asia Branch

#### Saikang Medical Technology Ltd

Sala Thammasop Subdistrict  
Thawi  
Watthana District  
Bangkok Thailand  
TEL: +66 642741046

### UK Branch

#### Uk Healthward International Limited

DEPT 111,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +31 (0) 2021111106

### Russia Branch

#### LLC 'SAIKANG'

454004, Russian Federation,  
Chelyabinsk,  
Universitetskaya Naberezhnaya,  
Russia  
TEL: +7 (919)300-00-06

### Mid East Branch


#### Saikang Medical Supplies (Mid East) LLC

Musaffah  
P.O Box 53200, Abu Dhabi,  
UAE  
Tel: +971 585710652

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

 +86 18021231901 ·  24/7

 export@sksurgitech.com

 www.sksurgitech.com



Company video



General catalogue