



# SKR-R06

## MONITOR CART

Registration code: SKM-A



SKSurgitech™  
by Saikang Medical



# SKR-R06

## FEATURES AND HIGHLIGHTS

### Product Introduction

The adjustable fixed base can perfectly match monitors of various sizes and is suitable for data monitoring in various clinical working environments.

### Technical Parameters

▪ Dimensions W×H	510×1000mm
▪ Storage frame size L×W×H	290×150×170mm
▪ Overall weight	7.3kg
▪ Caster design	3 inch medical grade silent wheel, pedal brake

### Optional Accessories

- Work surface



## Material

The table top and column are made of aviation-grade aluminum alloy, the base is made of one-time molding all-cast aluminum material, the whole car adopts high-voltage electrostatic spraying process, and imported environmentally friendly antibacterial plastic powder is used to ensure that the product meets the medical and health requirements, and the surface is matte and non-dazzling.



## Table Top

The table can be translated to both sides and has long screw fixing holes, which can be applied to monitors of various brands and sizes.



## Basket

With a height-adjustable storage basket for various measurement accessories required by the monitor.



## Column

The column is equipped with a manual lifting function, which is easy to rise and fall, and can be quickly adjusted to a suitable height and locked with a knob.



## Chassis

All-metal chassis, one-time cast aluminum molding, to ensure the center of gravity and level of the vehicle.



## Casters

3-inch medical-grade silent wheels, all with brakes, are flexible and durable after tens of thousands of kilometers of long-term durability tests.





## CONSTANT INNOVATION

Design

Research

Manufacturing

Maintenance

Service



## WORLDWIDE BRANCHES

### South America Branch

#### Saikang Medical Technology, Inc

W3400 Corrientes, Argentina  
TEL: +55 (11) 2758-0438  
Rua Universal, Hollywood,  
São Bernardo do Campo,  
Brazil  
TEL: +55 (11) 2758-0438

### East Africa Branch

#### Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,  
upper hill, nairobi, Kenya  
TEL: +254714084806

### South East Asia Branch

#### Saikang Medical Technology Ltd

Sala Thammasop Subdistrict  
Thawi  
Watthana District  
Bangkok Thailand  
TEL: +66 642741046

### UK Branch

#### Uk Healthward International Limited

DEPT 111,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +31 (0) 202111106

### Russia Branch

#### LLC 'SAIKANG'

454004, Russian Federation,  
Chelyabinsk,  
Universitetskaya Naberezhnaya,  
Russia  
TEL: +7 (919)300-00-06

### Mid East Branch



#### Saikang Medical Supplies (Mid East) LLC

Musaffah  
P.O Box 53200, Abu Dhabi,  
UAE  
Tel: +971 585710652

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

 +86 18021231901 ·  24/7

 export@sksurgitech.com

 www.sksurgitech.com



Company video



General catalogue