



# SKR-R07

## INFUSION PUMP CART

Registration code: SKM-A



SKSurgitech™  
by Saikang Medical

PRODUCT  
DATA SHEET



# SKR-R07

## FEATURES AND HIGHLIGHTS

### Product Introduction

The main structure of this trolley is made of high-quality aviation-grade aluminum alloy material, which is reliable in structure, flexible and convenient to use, and has the functions of an infusion pump and a syringe pump trolley. It is used in hospital wards, operating rooms, emergency centers, etc.

### Technical Parameters

▪ Dimensions L×W×H	580×490×1750mm
▪ Equipment mounting column diameter	25mm
▪ Storage frame size L×W×H	300×180×150mm
▪ Machine weight	13kg
▪ Caster design	3 inch medical grade silent wheel, pedal brake
▪ Rated voltage	220V
▪ Rated current	10A
▪ One for six power plug design	



## Material

The whole vehicle adopts medical grade aviation aluminum alloy, adopts high-voltage electrostatic spraying process, and uses imported environmentally friendly antibacterial plastic powder to ensure that the product meets the medical and health requirements, and the surface is matte and non-dazzling. The base is made of ABS medical grade plastic.



## I.V. Pole

The top of the column is equipped with a four-claw infusion stand, which can hang four different infusion bags at the same time.



## Equipment Column

The ultra-long one-piece equipment mounting column, with a diameter of 25mm, can perfectly match various infusion pumps and syringe pumps.



## Basket

There is a metal mesh basket on the column, which can store medicines and various groceries.



## Socket

The back is equipped with one-to-six power plugs, and is equipped with a power switch, which can support 6 devices at the same time.



## Casters

3-inch medical-grade silent wheels with brakes, after tens of thousands of kilometers of long-term durability tests, flexible and durable.





# CONSTANT INNOVATION

Design

Research

Manufacturing

Maintenance

Service



# WORLDWIDE BRANCHES

## South America Branch

### Saikang Medical Technology, Inc

W3400 Corrientes, Argentina  
TEL: +55 (11) 2758-0438  
Rua Universal, Hollywood,  
São Bernardo do Campo,  
Brazil  
TEL: +55 (11) 2758-0438

## East Africa Branch

### Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,  
upper hill, nairobi, Kenya  
TEL: +254714084806

## South East Asia Branch

### Saikang Medical Technology Ltd

Sala Thammasop Subdistrict  
Thawi  
Watthana District  
Bangkok Thailand  
TEL: +66 642741046

## UK Branch

### Uk Healthward International Limited

DEPT 111,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +31 (0) 2021111106

## Russia Branch

### LLC 'SAIKANG'

454004, Russian Federation,  
Chelyabinsk,  
Universitetskaya Naberezhnaya,  
Russia  
TEL: +7 (919)300-00-06

## Mid East Branch

### Saikang Medical Supplies (Mid East) LLC

Musaffah  
P.O Box 53200, Abu Dhabi,  
UAE  
Tel: +971 585710652

# SKSurgitech™

by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

+86 18021231901 · 24/7

export@sksurgitech.com

www.sksurgitech.com



Company video



General catalogue