



SKR-R09

INFUSION PUMP CART

Registration code: SKM-A



SKSurgitech™
by Saikang Medical

PRODUCT
DATA SHEET



SKR-R09

FEATURES AND HIGHLIGHTS

Product Introduction

The main structure of this trolley is made of high-quality aviation-grade aluminum alloy material, which is reliable in structure, flexible and convenient to use, and has the functions of an infusion pump and a syringe pump trolley. It is used in hospital wards, operating rooms, emergency centers, etc.

Technical Parameters

▪ Dimensions L×W×H	580×490×1750mm
▪ Equipment mounting column diameter	25mm
▪ Storage frame size L×W×H	300×180×150mm
▪ Machine weight	13kg
▪ Caster design	3 inch medical grade silent wheel, pedal brake
▪ Rated voltage	220V
▪ Rated current	10A
▪ One for six power plug design	



Material

The whole vehicle adopts medical grade aviation aluminum alloy, adopts high-voltage electrostatic spraying process, and uses imported environmentally friendly antibacterial plastic powder to ensure that the product meets the medical and health requirements, and the surface is matte and non-dazzling. The base is made of ABS medical grade plastic.



I.V. Pole

The top of the column is equipped with a four-claw infusion stand, which can hang four different infusion bags at the same time.



Equipment Column

The ultra-long one-piece equipment mounting column, with a diameter of 25mm, can perfectly match various infusion pumps and syringe pumps.



Basket

There is a metal mesh basket on the column, which can store medicines and various groceries.



Socket

The back is equipped with one-to-six power plugs, and is equipped with a power switch, which can support 6 devices at the same time.



Casters

3-inch medical-grade silent wheels with brakes, after tens of thousands of kilometers of long-term durability tests, flexible and durable.





CONSTANT INNOVATION

Design

Research

Manufacturing

Maintenance

Service



WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd

Sala Thammasop Subdistrict
Thawi
Watthana District
Bangkok Thailand
TEL: +66 642741046

UK Branch

Uk Healthward International Limited

DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +31 (0) 2021111106

Russia Branch

LLC 'SAIKANG'

454004, Russian Federation,
Chelyabinsk,
Universitetskaya Naberezhnaya,
Russia
TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC

Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

SKSurgitech™

by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

+86 18021231901 · 24/7

export@sksurgitech.com

www.sksurgitech.com



Company video



General catalogue