



# SKR-R10

## ECG MACHINE TROLLEY

Registration code: SKM-B



**SKSurgitech**<sup>™</sup>  
by Saikang Medical

Care you precisely

# SKR-R10

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ External size (LxWxH)	450x425x780mm
▪ Basket size	280x160x100mm
▪ Tray size	350x300mm
▪ Column diameter	50mm
▪ Base diameter	450mm



#### Design

The concise appearance design is suitable for the respiratory series, such as ECG, small monitoring equipment, etc.



#### Material

The column is made of aviation aluminum, the surface is anodized, it is not easy to rust and corrosion. Diameter 50mm. The base is made of aviation aluminum, the surface is painted, low center of gravity structure. Diameter 450mm. Load-bearing 100kg.



#### Basket

Metal immersion painting process, adjustable position.  
Size: 280\*160\*100mm.




#### Caster

5\*2.5 inch silent casters, flexible movement, 2 casters with brakes.



#### Handle

Integrated handle, streamlined design, smooth feel, ergonomic.



#### Tray

Aluminum tray with anodized surface treatment and firm installation. Size: 350\*300mm.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi,  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai,  
United Arab Emirates.  
Tel.: +20 128 8962497

**SKSurgitech™**  
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
☎️ +86 18021231901 · 🌐 24/7  
✉️ export@sksurgitech.com  
🌐 www.sksurgitech.com



Company video



General catalogue