



SKR054-AT

ANESTHESIA TROLLEY



SKSurgitech™
by Saikang Medical

Care you precisely

SKR054-AT

FEATURES AND HIGHLIGHTS

Technical parameters

- Size L850×W520×H950mm

Technical configuration

- Sliding Side Shelf 1pc
- Sharp Container 1pc
- Bust Baskets 2pc
- CPR Board 1pc
- Central Lock 1set
- File Box 1pc



Material

ABS body, metal ball-bearing rail.
Powder coated columns.



Table Top

Guardrails on top of trolley is made of stainless steel, table top with soft glass, avoid possibly falling down during movement.



Drawers

Inner with partitions can be well organized.



Storage Box

single row storage boxes.



Castor

High-quality medical silent universal castors, flexible and light to drive, smooth steering, wear-resistant, two with brake, can be easily braked in any state.



Drawer label card

The holder of the drawer is swallow tail design, outside with notice card slot, avoid liquid and dust access.



Accessories

Collapsible working table, file box;
2 dust basket, utility container;
Centralized lock;
Needle disposal holder.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🌐 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue