

SKR054-ET

EMERGENCY TROLLEY

Register code: SKM-A



SKR054-ET

FEATURES AND HIGHLIGHTS


Technical parameters

- External size (LxWxH) 850x520x960 mm
- Castors 125 mm



One-piece countertop

The one-piece countertop provides more comprehensive protection; the fully guarded handle makes it more comfortable to hold me; the raised grip shape makes the operation of the side-drawing board more comfortable.




Decorative gripper

Outer beveled decorative grip, better operating comfort and more three-dimensional shape.



Heavy duty construction

Wheels assembled on the same vertical line as four columns. More stable, even the drawer open inside full of medicine/tools.



Drawer

Drawer in size big & deeper!



Power supply

The power socket is movable and detachable, easy to use, located above the CPR board to convenient the use of oxygen cylinders



Inserted oxygen bottle holder

Inserted retractable oxygen bottle holder, belts fixed on trolley main body, space saving



Infusion stand

IV pole inserted into main body column.



Drawer holder

Card insert from downside, easy keep clean.



Defibrillator

The plate surface is made of virgin ABS material. The position of the I,V pole and the defibrillator can be exchanged for easy use.



Sorting trash can

Green+yellow sorted trash can.






CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™

by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
☎️ +86 18021231901 · 🕒 24/7
✉️ export@sksurgitech.com
🌐 www.sksurgitech.com



Company video



General catalogue