



SKR058-ET

EMERGENCY TROLLEY



SKSurgitech™
by Saikang Medical

Care you precisely

SKR058-ET

FEATURES AND HIGHLIGHTS

Material

The ABS material is an innovative engineering plastic. Antibacterial, moisture-proof, anti-corrosion, wear-resistant, non-toxic and tasteless, never rust. Washable, easy cleaning.



Design

The overall design and production are exquisite, the styling and color matching are reasonable, installation is simple, and the structure is stable and durable.



Drawer

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



Drawer compartment

There are transparent acrylic medicine grids inside the drawer, which can be used to classify and place medicines, and the medicine grids can be freely combined.

Technical parameters

Overall Dimensions	790W X 610D X 1000H mm
Defibrillator Stand	380×340mm
Work Board	340×250mm
Cpr Board	475 × 595 × 8.5mm
Packing Size	820X650x1150 mm
Product Weight	65 Kg
Caster	5"Two Casters With Brake, Two Without Brake
Work Board Size	340×250mm
Bracket Size	380×340mm


Technical configuration

- 20L Trash Can With Lid
- Oxygen Bottle Rack
- Adjustable Iv Stand
- Medicine Grid Horizontal And Vertical
- Slide Rail

Accessories


IV stand

Adjustable aluminum alloy IV stand.



Caster

At the bottom, it uses two 5-inch silent brake, anti-winding, high-end, silent casters, allowing for flexible movement and a high load-bearing capacity.




Drawer slide

The drawer uses high-strength three-section silent damping slides. When the drawer is closed to the final distance, hydraulic pressure is used to slow down the speed and reduce the impact force.




Working table

The main tabletop can have a replaceable soft pad that is dust-proof and easy to clean. The vehicle is equipped with a hidden sliding writing board that can be used for writing.



Sealed handle

The trolley is equipped with a sealed handle, which can lock five drawers at the same time.



Oxygen bottle rack

ABS oxygen bottle bracket on the left side.



Trash-can

Clamshell ABS trash can on the right side, Specification: 20L



CPR board

The ABS CPR board is installed on the back of the trolley body, which can be used during defibrillation rescue.

CPR board size: 475 * 595 * 8.5mm





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi,
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SKSurgitech™
by Saikang Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ export@sksurgitech.com

🌐 www.sksurgitech.com



Company video



General catalogue